PD VitalOs Cement®

Injectable and resorbable synthetic bone grafting cement
CalciphOs Technology

Ready-to-use sterile dual-paste cement conditioned in dual syringe
Single use product

Composition

Tricalcium Phosphate 36%, Monocalcium Phosphate 23%, Calcium Sulfate 11.5%, ultra pure water 27%, mineral excipient ad 100%.

Properties

PD VitalOs Cement® is a synthetic and biocompatible calcium phosphate cement intended for bone defect filling in dentistry. Once the cement paste has been injected, it can still be shaped and worked for a maximum of 2 minutes. It then takes approximately 10 minutes to set. When setting is completed, the strength of the cement is equivalent to that of the spongy bone. The final phase of the product is Brushite, a resorbable Calcium Phosphate. This cement promotes new autologous bone formation.

Indications

- > Bone regeneration in dental osseous defects.
- > Examples of indications: filling of periodontal infra-bony defects; sinus floor augmentation; filling of peri-implant defects.
- > For more information about the indications: www.vitalos.com

Directions for use

- > Prepare the site (bone defect) in the usual way and dry, avoiding as much as possible the presence of blood at the cement/bone interface (presence of blood between the cement and bone can hinder the anchoring of the cement to bone).
- > To prepare the syringe, remove the cap and replace it by the mixing tip.
- > Eject off-site the first cement drop (the mix is not optimal in this part of the product).
- > Inject PD VitalOs Cement® into the prepared site, shape to fit and complete the work as required. There is no need for a membrane.
- > Working time after injection of the paste: max. 2 minutes. In order to ensure a perfect setting and correct strength of the cement, do not work the paste longer than the recommended time.
- > Setting time to final hardness: approximately 10 minutes at 37°C.

Caution: after about 1 minute, the cement starts to set in the mixing tip and may block it. If more cement has to be injected (into the same site or on the same patient) simply place a new tip onto the syringe.

Note

PD VitalOs Cement® alone cannot provide stabilization of mobile teeth in severe periodontitis. A splinting treatment on such teeth is therefore indicated before injecting the cement (in direct contact with the tooth).

Packaging

Disposable sterile dual syringe + 2 mixing tips
Capacity of the syringe: 1 ml (~2 g)
1 plug to maintain humidity in the packaging (Note: humidity droplets may form on the surface of the packaging and of the product. This does not affect at all the quality of the product). Sterilized by gamma irradiation.

Storage

Store in its original protective packaging at room temperature between 5°C and 30°C, protected from direct light. If the protective packaging is damaged, discard the entire device as sterility can no longer be guaranteed

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Composition

Phosphate tricalcique 36%, phosphate monocalcique 23%, sulfate de calcium 11,5%, eau ultra-pure 27%, excipient minéral ad 100%.

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